



## **FIRST AID**

In case of contact with eyes, immediately flush eyes with large amounts of water for at least 15 minutes. In case of contact with skin, immediately wash skin with soap and flush with large amounts of water. If inhaled, remove to fresh air. If person not breathing, give artificial respiration. If breathing is difficult, give oxygen. If swallowed, wash out mouth with water provided person is conscious. Seek medical attention immediately. Wash all contaminated clothing before reuse.

## **-----REACTIVITY DATA-----**

**STABILITY:** Moisture-sensitive material.

**CONDITIONS TO AVOID:** Contact with ignition source or heat. Flammable liquid. Container may react with ignition source, heat, or fire conditions.

**INCOMPATIBILITY:** Strong oxidizing agents.

**HAZARDOUS DECOMPOSITION PRODUCTS:** Toxic fumes of carbon monoxide CO, carbon dioxide CO<sub>2</sub>, and nitrogen oxides.

**HAZARDOUS POLYMERIZATION:** None

## **-----ENVIRONMENTAL INFORMATION-----**

**SPILL OR LEAKAGE PROCEDURES:** Evacuate contaminated area. Shut off all sources of ignition. Wear self-contained breathing apparatus, rubber boots and heavy rubber gloves. Cover with dry lime or soda ash. Pick up and store in closed container and hold for appropriate waste disposal. Ventilate area well. Wash spill site after material pickup is complete.

**WASTE DISPOSAL:** Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber. Observe and consult state, local, or federal EPA regulations for proper disposal and environmental regulations.

## **-----PROTECTION AND PRECAUTIONS-----**

**VENTILATION REQUIREMENTS:** Mechanical exhaust required.

**RESPIRATORY PROTECTION:** NIOSH/MSHA approved respirator.

**PROTECTIVE GLOVES:** Chemical resistant and safety gloves.

**EYE/FACE PROTECTION:** ANSI approved safety goggles

**HANDLING AND STORAGE:** Keep container tightly closed. Keep away from heat, sparks, and open flame. Store in a cool dry place. Do not breathe vapor. Irritant.

**OTHER PRECAUTIONS:** Flammable. Irritant. Irritating to eyes, respiratory systems and skin. Keep away from source of ignition. No smoking. Keep container tightly closed in a cool well ventilated place. Wear suitable protective clothing. Moisture sensitive. Store under nitrogen.

-----  
The information herein is believed to be accurate and reliable as of the date compiled. However, Chemat Technology, Inc. makes no representation, warranty or guarantee of any kind with respect to the information on this data sheet or any use of the product based upon this information.

DATE PREPARED: April 2004